

MiZuki

MiZuki



Micro Catheter




KANEKA
MEDICAL
PRODUCTS



Guidewire Support Catheter for real world PCI

Product Overview



Features

1 Advanced Hydrophilic Coating

Our new technology in the hydrophilic coating realizes excellent trackability and crossability for the success in the procedures. (in comparison with data of our existing products)



2 Product lineup with 2 types

Two types with different tip stiffness enlarge the product applicability to various cases.

3 Advanced shaft design

Fluoro-resin inner surface for lubricity reinforced with braiding and the smooth transition of the shaft stiffness achieves guidewire maneuverability, trackability and crossability.

Product Specifications

Product Name	Catalogue Number	Catheter Effective Length(cm)	Shaft Outer Diameter(Fr/mm)			Shaft Inner Diameter (inch/mm)	
			Tip	Distal	Proximal	Distal	Proximal
 (Standard)	MZ14-135BM	135	1.8/0.60	1.8/0.60	2.5/0.84	0.018/0.45	0.022/0.55
	MZ14-150BM	150					
 (FX type)	MZ14-135TM	135	1.7/0.58	1.7/0.58		0.017/0.42	
	MZ14-150TM	150					



Do not re-use



Sterilized using ethylene oxide



Consult instructions for use

Manufacturer

KANEKA CORPORATION

3-18, 2-Chome, Nakanoshima, Kita-ku, Osaka-city, OSAKA, 530-8288 JAPAN

TEL. +81-6-6226-5256 **CONTACT** TEL. +81-3-5574-8136