


---


Directional  
Coronary  
Atherectomy

---

Directional Coronary Atherectomy  
**ATHEROCUT**™



## Catheter design



Model No.	Size	Recommended Vessel Diameter	Compatible Guiding Catheter	Compatible Guidewire
NDCA-S9	S	3.0-3.4mm	7F (inner diameter $\geq$ 2.0mm)	0.014inch (0.36mm)
NDCA-M9	M	3.5-3.9mm		
NDCA-L9	L	4.0-4.4mm		

### MDU (Motor Drive Unit)



Model No.  
NDCA-MDU


### Inflation device



Model No.  
NDCA-Indefrator

## Compliance chart

Inflation Pressure (atm)	Size: S (for vessel diameter 3.0-3.4mm)	Size: M (for vessel diameter 3.5-3.9mm)	Size: L (for vessel diameter 4.0-4.4mm)
1	3.09	3.45	3.95
2	3.31	3.71	4.21
3	3.45	3.86	4.39
4	3.56	3.98	4.52



VASCULAR DIV.

[http://med.nipro.co.jp/en\\_products](http://med.nipro.co.jp/en_products)



NIPRO CORPORATION  
3-9-3, Honjo-Nishi, Kita-Ku, Osaka 531-8510, Japan

- ▶ Please consult your sales representative regarding stock availability.
- ▶ Please refer to the Instructions For Use for further details regarding product use and specifications.

[Manufacturing facility]  
NIPRO CORPORATION ODATE FACTORY  
8-7, Hanukiyachi, Niida, Odate-shi, Akita 018-5794 Japan

DCA-1810-E1