




Ergonomic design, Reliable performance

Complying with the most stringent of test parameters, the GM-30 is representative of the emphasis NIPRO places on producing reliable, high quality products, while the ergonomic design ensures ease-of-use.

GM-30

Inflation Device

Features



Pressure gauge

Maximum pressure reading:
3.0MPa (30atm / 441psi)

Fluorescent display allows for viewing
in darkened conditions



Release lever

Ergonomic design that fits comfortably
in the hand to allow quick release with
minimum effort



Handle

Designed for easy rotation and handle
extraction

Product Information

Maximum Pressure Reading	3.0MPa (30 atm / 441 psi)
Maximum Volume Dispensed	25 mL
Accuracy	±3 %FS

Ordering Information

Model Number	GM-30F
--------------	--------

► Please consult your sales representative regarding stock availability.